

#### 7.1.10

##### meeting minutes

run controls on digestions, check emails for email from francis about what controls to run

Run digestions at 37 degrees for 15 min, one for an hour, one for overnight to see which work best. Most protocols say 1 hour for digestions and ligations. Reaction mixes are fine.

One with no enzyme, one with no DNA, ...steve is sending an email with all this info

We should aliquot out all our buffers, dNTPs.

Whenever we get ladders, make sure we dilute, aliquot . can use 1 ul of ladder for each gel for clear bands.

Start saving electroporation cuvettes, we can re-use them.

Can register computer to printer on north side. Go to IGB website. The IP address you need is on the printer behind the secretary's desk.