09. (September) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Theresa Keil WEDNESDAY, 11/9/2019

• prepared RIPA buffer 2x after paper c-Src

150ml 2x RIPA

20mM Tris-Hcl -> 0.479g

300mM NaCl -> 2,63g

2mM EDTA -> 0.088g

0.2% SDS

2% Triton

2% Sodiumdeoxycholate

(1:50 Protease-Inhibitor-Cocktail is added freshly)

Table4												
	Α	В	С	D	E	F	G	Н	I	J	K	L
1	Marker	Exo His MCP SN	Exo MCP SN	Exo L7Ae His SN	Exo L7Ae	VLP Biotin MCP SN		Marker	VLP Biotin MCP SN	VLP SN His	VLP SN Mock	Exo His MCP SN

file:///tmp/tmpGJR1v8.html