

09. (September) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Anja

MONDAY, 2/9/2019

Joshi

Cell culture: HEK293T splitting and seeding (2x 6-well plate)

- 750.000cells/well, 2.5ml full DMEM/well
- splitted 400ul

Joshi & Annika

Restriction Digest (according to Restriction Digest)

Ligation (as previously described)

Transformation (according to Chemical transformation)

- E16, 17 for Glycerol stocks

Anja

Cell culture: HEK293T splitting and seeding (2x 96-well plate)

Origin cells, 25.000 cells/well, 150 ul full DMEM/well