

08. (August) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Johanna Wallner

THURSDAY, 15/8/2019

Alejandro

cell culture: Medium exchange 96-well plate HiBit assay + 24-well plate qPCR and Purification

- 100 μ L medium from the 96-well plate (HiBit assay) were exchanged (11:30 \pm 48 h post transfection)
- 600 μ L medium from the 24-well plate with 8 wells (qPCR + Affinity purification samples) were exchanged (11:30 \pm 49.5 h post transfection)