

06. (June) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Johanna Wallner

SUNDAY, 2/6/2019

Cell culture medium preparation:

Medium was mixed according to cell culture medium preparation protocol.

(500 ml Advanced DMEM (Gibco) + 50ml FBS (10%) + 6 ml GlutaMAX (Gibco) + 6 ml Pen-Strep (Gibco)).


 Cell Culture Medium Preparation

Cell splitting:

HEK 293T cells were splitted according to cell splitting protocol.

Measurements from cell counting : 3.4×10^6 cells/mL.

Passage 50.

 Cell Splitting