

07. (July) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Theresa Keil

WEDNESDAY, 31/7/2019

continued with BS from 30.07 VLPs and EXO stripped after protocol

made two new aliquots of secondary ab goat anti mouse HRP

each: 1.5000 in 5 % milk in 1xTBS

used old aliquot for exo-stripped and one new for VLPs 2nd one frozen stored at -20 °C

threw away sec AB from 26.06

- finished like WB protocol
- continued with VLP blot from 30.07.19
- stripping
- used primary ab: p24: 10 µl in 5 ml 5% milk TBS

continued RNA isolation from cells and exosomes from 30.07

normal continued protocol according to Trizol user guide stopped after resuspending in 70 % ethanol

stored at -20 °C