## 10/6/2019

## 04. (April) 2019

Project: iGEM\_Munich2019 Shared Project

Authors: Sarah Brajkovic

TUESDAY, 30/4/2019

**Transfection**: Fussenegger Plasmids

Cells were transfected in a 12-well plate at 70%-90% confluency as shown below.

Wells were transfected with Lipofectamine according to the lipofectamin 3000 transfection protocol.

300419_12-well plate for transfection				
	1	2	3	4
А	pSA462	pSA462	Booster	Booster
В	pSA465	pSA465	Booster pSA462 pSA465	Booster pSA462 pSA465
С	pSA462 pSA465	pSA462 pSA465	Ctrl: PBS	Ctrl: PBS

file://tmp/tmpi2r723.html