

06. (June) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Anja

TUESDAY, 11/6/2019

Transfection (10cm dish)

Table3			
	Component	Amount	Volume
1	DNA (E6*)	7,5ug	6,88ul
2	DNA (E12*)	7,5ug	21,25ul
3	PEI		45ul
4	Optimem		1000ul

Incubation 15min at RT

*for concentrations see 4.6.2019

MaxiPrep overnight inoculation preparation (following the protocol)