

09. (September) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Anja

MONDAY, 16/9/2019

Anja

Cell culture: HEK293T splitting and seeding

- 96-well plate, 25.000 cells/well, 150 ul full DMEM/well (EXOs & VLPs, for CTX-general & CTX-killing assay in 3 time conditions - 24, 48, 72h, 6 plates)

Alejandro

EXOs HiBiT Assay for HEK293T(see protocol)

Joshi

Cell culture: HEK293T splitting and seeding

- 2x 6-well plates seeded
- 2x p75 flasks splitted (EXOs and VLPs)