09. (September) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Anja

TUESDAY, 10/9/2019

Joshi

Cell culture: Transfection of HEK293T

| 10092019_Transfection_E | | | | | |
|-------------------------|------------|------------|------------|--|--|
| | 1 | 2 | 3 | | |
| Α | AL, MeL | AM, MeM | AL, MeM | | |
| В | AM, MeL | М | 1 | | |

| 10092019_Transfection_E | | | | | |
|-------------------------|-------------------|-------------------|---|--|--|
| | 1 | 2 | 3 | | |
| Α | CC, AL, MeL | CC, AM, MeM | 1 | | |
| В | 1 | 1 | 1 | | |

Preparation for sequencing: E18

Anja

MidiPrep of E18 (see protocol)

• eluted in 300ul

Concentration:

-E18 = 405.2ng/ul (prepared 100ng/ul dilution)

Moritz

Cell culture: HEK293T splitting and seeding

One 96-well plate was seeded to try the cytotoxicity killing assay with Zeocin/PBS diluted DMEM (28.000 cells/well, 150ul).

file:///tmp/tmpRXagDW.html