

07. (July) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Johanna Wallner

SATURDAY, 6/7/2019

HiBit assay:

- The following samples were thawed from the 5th of July, then diluted 1:10 in PBS and the HiBit Assay was repeated (50µL sample + 50µL reaction mix; Extracellular detection system).
 - 72h, cell lysate, a-c
 - 72h, lysed VLPs, a-c
 - 72h, unlysed VLPs, a-c
- Incubation on shaker : 5min
- FinalTriton X-100 concentration : 0.05%.