## 09. (September) 2019

Project: iGEM\_Munich2019 Shared Project

Authors: Johanna Wallner

MONDAY, 23/9/2019

## <u>Johanna</u>

## SDS-PAGE Purification

- repeat the SDS-PAGE from saturday
- dilute SN and FT 1:10
- V8 (SN, FT, W, E), MCP (SN, FT, W, E), L7Ae (SN, FT, W, E)



## Cell culture

- seeding:
  - 24 wells of a 96-well plate were seeded: 25000 cells / well in 150 μL medium
  - o a full 6-well plate was seeded: 750000 cells/well in 2.5 mL medium
  - o time: 11:30
- Splitting: passage 37

file:///tmp/tmp5p2acE.html