06. (June) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Theresa Keil
THURSDAY, 6/6/2019

Exosome precipitation

samples were taken out of the fridge after 24 h centrifugation at 2,000xg for 30 min at 4 °C, pour off supernatant 8,000 PEG taken off supernatant 35,000 PEG resuspension in PBS 100 µI negative sample, 500 µI for sample

His-Tag-Purification of Exosomes using Agarose Ni-NTa resin.

file:///tmp/tmpnTZHKV.html