

08. (August) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Anja

THURSDAY, 22/8/2019

Johanna, Sarah

His-Tag Purification (Protocol: His-Tag-Purification (Gravity Flow, 2 mL sample)) (repeat, see 18.08.2019 entry)

Constructs Tested:

- E11: CD63-MCP
- E17: CD63_Ser161_6xHis-MCP
- E12: Lamp2b_6xHis
- V30: mock

with adjustments:

- Wash Buffer: 15mM imidazole instead of 25mM