

# 06. (June) 2019

---

**Project:** iGEM\_Munich2019 Shared Project

**Authors:** Johanna Wallner

FRIDAY, 7/6/2019

---

## Cell splitting:

HEK 293T cells were splitted according to cell splitting protocol.

Measurements from cell counting :  $1.8 \times 10^6$  cells/mL.

Passage 52.

 Cell Splitting

## Cryostocks of V1, V4, and V5

750 µl of bacterial overnight cultures were mixed with 750 µl autoclaved glycerol in Cryotubes and frozen at  $-80^{\circ}\text{C}$ .