

07. (July) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Sarah Brajkovic

FRIDAY, 12/7/2019

Anja

HEK293T splitting_(p16)

Cells were splitted (500ul into a new p75 flask)

MaxiPrep of V15/E9 (Promega PureYield Maxi kit)

Concentration : average 103.3ng/ul in approx 1.2ml, ratios > 1.5, sent for sequencing

Joshi

Fluc-Reuptake Assay (see protocol)

Cells were analysed for fluc content 24h after media exchange. Separate samples were taken for PCR (fluc).