

Modified Compotent Cells Procedure

Introduction

Get started by giving your protocol a name and editing this introduction.

Materials

- › Dh5a cells
- › Microfuge tubes
- › ice-water bath
 - › Styrofoam bucket
 - › float(s)
- ›

Procedure

Preliminary Steps

- ✓ 1. Start liquid culture of Dh5a cells (single colony) from plate.
Grow overnight
- ✓ 2. Get _ pre-chilled MCF tubes from fridge.
- ✓ 3. Prepare ice-water bucket.

Main Procedure

- ✓ 4. Evenly transfer Dh5a cells into MCF tubes.
- ✓ 5. Centrifuge for 10 minutes at 6000 rpm
Use the smaller clinical centrifuge