

AMMONIUM SULPHATE PRECIPITATION

15,05 g ammonium sulphate were added to 50 mL supernatant (L1-sAP-RBD-Strep8His) -> 50% ammonium sulphate saturation

The mixture was stirred on ice for 20 min.

https://www.sheffield.ac.uk/polopoly_fs/1.711327!/file/Protocols_and_tips_in_protein_purification.pdf

<https://www.encorbio.com/protocols/AM-SO4.htm>

<https://www.ucl.ac.uk/~ucbcdab/enzpur/amso4.htm>