
DATASHEET: THEOPHYLLINE RIBOSWITCH 1

Part Number	Part Name	Construct	Length
BBa_K537003	Theophylline riboswitch 1-Venus	throughlylling illow olds h 1 (this) lyvenus (6.5 / 2.035)	780 bp
BBa_K537009	Promoter-Theophylline riboswitch 1-Venus- Double terminator	Storag promote	9601 bp
BBa_K537001	Theophylline riboswitch 1-CheZ N-fusion	Theopylline dopwitch 1 (nRS1)-Cne2 (NS3)001)	698 bp
BBa_K537011	Promoter-Theophylline riboswitch 1-CheZ Venus-Double terminator		1583 bp

Theophylline riboswitch 1: Summary

Using a high-throughput screen(Lynch et al., 2007) and a high-throughput selection (Topp and Gallivan, 2008) in *E.coli*, the Gallivan lab has been able to develop a powerful theophylline-dependent riboswitch. It is worth noting that the motility-based screen (Topp and Gallivan, 2008) made use of a theophylline riboswitch library activating the expression of CheZ, the motility factor in *E.coli* which allows for running (as opposed to tumbling) movement. Therefore, activation of CheZ allows bacteria to move along the theophylline gradient. The best riboswitch identified in these earlier screens was the theophylline riboswitch clone 8.1 (Figure 1) which displays 36-fold expression activation ratio.

8.1 Theophylline RS (Topp&Gallivan JACS 2007)






Figure 2: The fluorescence produced by Riboswitch 1 CheZ Venus (BBa_K537011) (A) and Riboswitch 1 Venus (BBa_K537009) (B) at different theophylline concentrations over time. Bacterial cell cultures were transformed with the respective constructs and grown to the mid-log phase of growth from a seeding culture. The cultures were excited at 514 nm and emission intensity was detected at 528 nm using a Jasco FP-6300 spectrofluorometer. Activation of the riboswitch can be seen in the presence of theophylline. The presence of the CheZ gene may confer some structural stability that enhances the activation when compared to the construct with venus alone. The activation seen may be sufficient to restore motility in CheZ deficient *E. coli* cells. Data is represented in triplicate ± SEM. * p < 0.05.

Fluorescence microscopy




Figure 3: *E.coli* CheZ mutants which were transformed with the Riboswitch 1 Venus (BBa_K537009) and Riboswitch 1 CheZ Venus (BBa_K537011) constructs. The bright field (BF) images in the left column depict all bacterial cells in the field. The Venus fluorescence images in the right column depict bacterial cells which emitted fluorescence. In the absence of theophylline, almost no fluorescence occurred (as with parental controls - left panels). Upon the addition of theophylline at a concentration of 1.5mM, many of the cells emitted fluorescence showing activation of the theophylline riboswitch 1.

Motility Assay




Figure 4. The distance travelled by $E.\ coli$ containing either Riboswitch 1 CheZ Venus (BBa_K537011) or Riboswitch 1 Venus (BBa_K537009) (C) at different theophylline concentrations over time. The observed trend in A can be seen on semi-solid agar plates containing 1.0mM theophylline with cells transformed with Riboswitch 1 CheZ Venus (B) and Riboswitch 1 Venus (C), where the halo surrounding the point of inoculum in B is indicative of regained motility. Data is represented in triplicate \pm SEM. *p < 0.05.




Figure 5. A graph to show the chemotaxis index of $\Delta \text{CheZ } E.\ coli$ transformed with Riboswitch 1 CheZ Venus, Riboswitch 1 Venus and RBS CheZ Venus (BBa_K537011, BBa_K537009, and BBa_K537013). The *Chemotaxis Index* of 100 for BBa_K537009 indicates bacteria travelled equally towards 0mM theophylline as 2.0mM theophylline, suggesting no attraction towards theophylline. Conversely, the chemotaxis index of 230 for BBa_K537011 increased movement towards 2.0mM theophylline compared with 0mM theophylline, strongly indicating riboswitch-controlled chemo-attraction. Data is represented in triplicate \pm SEM. * p < 0.05.

References

Lynch S.A. Desai S.K., Sajja, H.K., and Gallivan J.P.A high-trhoughput screen for synthetic riboswitches reveals mechanistic insight into their function. 2007, ChemBiol 14:173-18

Topp S. and Gallivan J.P. Random walks to synthetic riboswitches – a high throughput selection based on cell motility. 2008, ChemBiochem 9:210-213