

7/25

- Started making competent cells with DH5α.
- loaded gel (digests)

Lane 1 (top): [ladder] K081012 #1 | K081012 #2 | K081012 #3 | J23100 #2

Lane 2 (bottom): [ladder] J23102 #2 | J23104 #1 | LB digest #1 | LB digest #2 | I721002 #1 | I721002 #2 | I721002 #3
low DNA

- Transformation of ligations: K174015 + kan backbone using new and old competent cells (8 samples total) K824008 A + Chlor K824008 B + Chbr K824012 + Chlor

* note, transformations sat overnight in LB and SOC and were plated in the morning (full ~~1 ml~~ 1 ml without draining supernatant)