

PHOSPHATE SENSOR MODULE: LAB NOTEBOOK WEEKLY SUMMARY

Week 1 (1/6/2015-5/6/2015)

1. Summer training, mainly for recalling the basic techniques
2. Discussion for the plan of the main project

Week 2 (8/6/2015-12/6/2015)

1. Transformation of pSB1C3-BBa_I0500, pSB1C3-BBa_B0015, pSB1C3-BBa_B0030, pSB1C3-BBa_K381001 and pSB1C3-BBa_E0840
2. Inoculation of pSB1C3-BBa_I0500, pSB1C3-BBa_B0015, pSB1C3-BBa_B0030, pSB1C3-BBa_K381001 and pSB1C3-BBa_E0840
3. Plasmid extraction of pSB1C3-BBa_I0500, pSB1C3-BBa_B0015, pSB1C3-BBa_B0030, pSB1C3-BBa_K381001 and pSB1C3-BBa_E0840
4. Phosphorylation and annealing of *P_{phoA}*
5. Digestion of pSB1C3-BBa_J04450
6. Ligation of *P_{phoA}* with pSB1C3 and pSB1C3-BBa_I0500 with pSB1C3-BBa_B0030

Week 3 (15/6/2015-19/6/2015)

1. Transformation of pSB1A2-BBa_B0030, pSB1A2-BBa_E0840, pSB1AK3-BBa_B0015
2. Phosphorylation and annealing of *P_{phoA}*
3. Digestion of pSB1C3-BBa_J04450

Week 4 (22/6/2015-26/6/2015)

1. Colony PCR for checking the identity of pSB1C3-*P_{phoA}* candidate 1 and 3 by VF₂ and VR primers
2. Restriction check of pSB1C3-*P_{phoA}* candidate 1 and 3 using *PvuII*

Week 5 (29/6/2015-3/7/2015)

1. Digestion of pSB1C3-BBa_B0030.
2. Digestion of pSB1C3-BBa_J04450.
3. Overnight ligation of *P_{phoA}* with pSB1C3 at 16°C.
4. Transformation of pSB1C3-*P_{phoA}*.

5. Colony PCR for checking the identity pSB1C3-*P_{phoA}* candidate colonies. Identified candidate 1-17 of pSB1C3-*P_{phoA}*.
6. Digestion of pSB1C3-BBa_B0030
7. New plan for *P_{phoA}*, instead of having the construct of pSB1C3-*P_{phoA}* for sequencing, rather, ligate *P_{phoA}* with pSB1A2-BBa_E0840 for sequencing.
8. Digestion of pSB1A2-BBa_E0840
9. Overnight ligation of pSB1A2-BBa_E0840 with *P_{phoA}* at 16°C.

Week 6 (6/7/2015-10/7/2015)

1. Transformation of pSB1A2-*P_{phoA}*-BBa_E0840.
2. Identified pSB1A2-*P_{phoA}*-BBa_E0840 candidates by restriction check using *PvuII*.
3. Miniprep of pSB1A2-*P_{phoA}*-BBa_E0840 candidate 1-5.
4. Ligation of pSB1A2-BBa_E0840 with *P_{phoA}*.
5. Send pSB1A2-*P_{phoA}*-BBa_E0840 candidate 5 for sequencing.
6. Inoculation and streak pSB1C3-BBa_B0030, pSB1AK3-BBa_B0015

Week 7 (13/7/2015-17/7/2015)

1. Send pSB1A2-*P_{phoA}*-BBa_E0840 for sequencing.
2. Confirmed pSB1A2-*P_{phoA}*-BBa_E0840 from the sequencing result.
3. From sequencing result of pSB1C3-BBa_B0030, it showed extra nucleotides before and after the BBa_B0030 sequence, which is inconsistent with the sequence from the iGEM Parts Registry.
4. Transformation of pSB1C3-BBa_B0032.
5. Miniprep of pSB1C3-BBa_B0032 and pSB1C3-BBa_B0030.
6. Digestion of pSB1C3-BBa_B0032.
7. Glycerol stock for pSB1AK3-BBa_B0015.
8. Functional assay for *P_{phoA}*.
9. Characterization of *P_{phoA}*.
10. Prepared modified M9 minimal medium for *P_{phoA}* characterization with phosphate being replaced.

Week 8 (20/7/2015-24/7/2015)

1. Inoculate pSB1A2-*P_{phoA}*-BBa_E0840 in M9 minimal medium with 0μM, 10μM, 500μM and

1000 μ M of PO₄ (Using Tris hydrochloride to replace the original PO₄ in M9)

2. Check OD and GFP expression

3. Digestion and gel purification of pSB1A2-BBa_E0840 (for ligating with pSB1A2-BBa_J23101 to act as a positive control)

Week 9 (27/7/2015-31/7/2015)

1. Characterization of pSB1A2-*P_{phoA}*-BBa_E0840 in M9 minimal medium with 0 μ M, 30 μ M, 50 μ M, 100 μ M, 150 μ M, 200 μ M, 250 μ M PO₄³⁻.

2. Digestion of pSB1A2-BBa_J23101-J61002 (constitutive promoter-RFP reporter)

3. Ligation of BBa_J23101 (constitutive promoter) and BBa_E0840 (GFP reporter)

4. Streak plates for pSB1A2-*P_{phoA}*-BBa_E0840 candidate 5, pSB1A2-E0840 and pSB1A2-BBa_J23101-E0840

Week 10 (3/8/2015-7/8/2015)

1. Characterization of pSB1A2-BBa_*P_{phoA}*-E0840.

Week 11 (10/8/2015- 14/8/2015)

1. Passed to Rice University.

Week 12 (17/8/2015- 23/8/2015)

1. Continued the work of phosphate sensor from Rice University.

6. Characterization on *P_{phoA}*.

7. Characterization on *P_{phoBR}*.

8. Characterization on *P_{amn}* promoter.

Week 13 (24/8/2015-30/8/2015)

1. Characterization on *P_{phoA}* promoter.

2. Characterization on *P_{phoBR}* promoter.

3. Characterization on *P_{amn}* promoter.

Week 14 (31/8/2015-6/9/2015)

1. Characterization on *P_{phoA}* promoter.

2. Characterization on P_{phoBR} promoter.
3. Characterization on P_{amn} promoter.