

08. (August) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Johanna Wallner

MONDAY, 12/8/2019

Alejandro

Cell culture

- splitting
 - cell concentration 5 mL suspension: 3.1×10^6 cells/mL
 - new motherplate: passage 25, 1 Mio cells
- seeding
 - 11 wells of 24-well plates were seeded: 150000 cells in 750 μ L for Purification and qPCR
 - 48 wells of a 96-well plate were seeded: 25000 μ L in 150 μ L for HiBit
 - Time: 12:00 a.m.