

06. (June) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Johanna Wallner

TUESDAY, 25/6/2019

VLP harvesting_:

- harvesting was done after 51h and not 48h
- the culture medium was completely yellow. Seeding at 30.000 cells per well was too much. Confluency was at 100-110%, but there were no dead cells visible
- the harvesting was carried out according to the 150µL culture medium procedure. All supernatants were tested and the cytosolic content of P2A-/-, CAP,2:1,CAR and V8.
- the cell lysis was done for 45min instead of 30min
- samples were frozen at -20°C before the incubation at 60°C step