

08. (August) 2019

Project: iGEM_Munich2019 Shared Project

Authors: Johanna Wallner

WEDNESDAY, 21/8/2019

Alejandro

cell culture:

- medium exchange from the 3 6-well plates
- 1.6 mL were taken out (650 µL still remain) and 2.6 mL new medium were added --> total volume: 3.2 mL
- finished at 16 :30, 27 h post transfection

transformation of the new plasmid

- BioBrick BBa_K2170002 (Munich iGEM Team 2016) was obtained from Volker Morat
- 20 ng plasmid were used to transform NEB stable E. coli with the standard protocol
- 80 µL (x2) bacterial suspension was plated on an Ampicillin plate and a chloramphenicol plate (25 µg/ml)